

Panbio™ COVID-19 Ag Rapid Test manufactured by **Abbott** is available for procurement through UNICEF Supply Division. The test is available in 2 configurations, both offered at the same FCA-based price US \$2.25 per test.

Technical and commercial description of the product is presented in Table 1 below.

Table 1: Panbio™ COVID-19 Ag Rapid Test

UNICEF Material Number	Product Description	Product Code	Pack Size	Price per Pack, USD	FCA Lead Time, <i>indicative</i>	FCA Delivery Port
S0003039	Panbio COVID-19 Ag Rapid Test, kit/25T Kit content: 25 x Test devices with desiccant in individual foil pouch. 1 x Buffer bottle (9 mL). 25 x Sterilized nasopharyngeal swabs for sample collection. 25 x Extraction tubes. 25 x Extraction tube caps. 1 x Positive control swab. 1 x Negative control swab. 1 x Tube rack. 1 x Quick reference guide (Nasopharyngeal). 1 x IFU.	41FK10	25 tests	56.25	3 weeks	FCA ADK WH, South Korea, or FCA DHL Nijmegen, the NL.
S0003040	Panbio™ COVID-19 Ag Rapid Test Device (NASAL) kit/25 tests. <i>(No 2D barcode)</i> Kit content: 25 x Test devices with desiccant in individual foil pouch. 1 x Buffer bottle. (9 mL) 25 x Sterilized nasal swabs for sample collection. 25 x Extraction tubes. 25 x Extraction tube caps. 1 x Positive control swab. 1 x Negative control swab. 1 x Tube rack. 1 x Quick Reference Guide. 1 x IFU.	41FK11	25 tests	56.25	3 weeks	Same as above

Product Reference Note:
Panbio™ COVID-19 Ag Rapid Test by Abbott

Intended use:

Panbio™ COVID-19 Ag Rapid Test Device is an in vitro diagnostic rapid test for the qualitative detection of SARS-CoV-2 antigen (Ag) in either human nasopharyngeal or nasal swab specimens from individuals who meet COVID-19 clinical and / or epidemiological criteria. Panbio™ COVID-19 Ag Rapid Test Device is for professional use only and is intended to be used as an aid in the diagnosis of SARS-CoV-2 infection.

- The test provides preliminary test results and they cannot be used as the sole basis for treatment or other management decisions.
- The test is not intended to be used as a donor screening test for SARS-CoV-2.

The **Instructions for Use** for the test can be found in the following link: [IFU](#)

WHO Interim Guidance on the Use of COVID-19 Antigen-detecting Tests can be found on the following link: [WHO interim guidance](#)

Storage and stability conditions:

Store at 2-30 °C.

Shelf life:

24 months

Real-time Data Collection and Visualization:

Panbio™ COVID-19 Ag Rapid Test Device is compatible with Sympheos™ System.

Should you have any questions regards to ordering please contact Tatsiana Ptashnik tptashnik@unicef.org

Should you have any technical questions regards to the test please contact Wandani Sebonego wsebonego@unicef.org

Product Reference Note:

Panbio™ COVID-19 Ag Rapid Test by Abbott